

Lacrimal Cannula (19G) (9977D)

Images



Specifications

Property	Value
Overall Length (mm)	41
Gauge (mm)	19G
Luer Lock (mm)	EN ISO 20594

Other Details

- This product can be viewed online at http://www.eyeline.org/products/view/lacrimal_19g.

Lacrimal Cannula (26G) (9978D)

Images



Specifications


Property	Value
Overall Length (mm)	41
Gauge (mm)	26G
Luer Lock (mm)	EN ISO 20594

Other Details

- This product can be viewed online at http://www.eyeline.org/products/view/lacrimal_26g.

Sub-Tenons Cannula (9979D)

Images



Specifications


Property	Value
Overall Length (mm)	37.5
Luer Lock (mm)	EN ISO 20594

Other Details

- This product can be viewed online at <http://www.eyeline.org/products/view/sub-tenons>.

Angled/Rycroft Hydrodissection Cannula (9980D)

Images



Specifications

Property	Value
Overall Length (mm)	30
Gauge (mm)	25G
Luer Lock (mm)	EN ISO 20594

Other Details

- This product can be viewed online at http://www.eyeline.org/products/view/hydrodissection_rycroft.

Helsinki Hydrodissection Cannula (9981D)

Images



Specifications

Property	Value
Overall Length (mm)	35
Gauge (mm)	26G
Luer Lock (mm)	EN ISO 20594

Other Details

- This product can be viewed online at http://www.eyeline.org/products/view/hydrodissection_helsinki.

Bolger Hydrodissection Cannula (9982D)

Images



Specifications


Property	Value
Overall Length (mm)	33
Gauge (mm)	26G
Luer Lock (mm)	EN ISO 20594

Other Details

- This product can be viewed online at http://www.eyeline.org/products/view/hysrodissection_bolger.

Angled Cystiotome (9983D)

Images



Specifications


Property	Value
Overall Length (mm)	23
Gauge (mm)	27G
Luer Lock (mm)	EN ISO 20594

Other Details

- This product can be viewed online at http://www.eyeline.org/products/view/cystiotome_straight.

Curved Cystiotome (9984D)

Images



Specifications


Property	Value
Overall Length (mm)	25
Gauge (mm)	25G
Luer Lock (mm)	EN ISO 20594

Other Details

- This product can be viewed online at http://www.eyeline.org/products/view/cystiotome_curved.

Simcoe Irrigation/Aspiration Cannula (9985D)

Images



Description

Dual-bore irrigation/aspiration cannula & luer syringe fitting with 250mm silicone tubing.

Specifications

Property	Value
Overall Length (mm)	32
Luer Lock (mm)	EN ISO 20594

Other Details

- This product can be viewed online at http://www.eyeline.org/products/view/simcoe_cannula.