

Chalazion Lid Clamp (9970D)

Images



Specifications


Property	Value
Material	Stainless Steel
Overall Length (mm)	91
Max Handle Width (mm)	7.5
Clamp Ring Outer Diameter (mm)	10.5
Clamp Ring Inner Diameter (mm)	8.2
Max Clamp Opening (mm)	9

Other Details

- This product can be viewed online at http://www.eyeline.org/products/view/chalzion_clamp.

Chalazion Double-Ended Curette (9972D)

Images



Specifications

Property	Value
Material	Stainless Steel
Overall Length (mm)	135
Handle Style	Square Section, Satin
Smaller Spoon Diameter (mm)	1.7
Larger Spoon Diameter (mm)	2.6

Other Details

- This product can be viewed online at <http://www.eyeline.org/products/view/curettes>.

Scalpel No3 Handle & Blade No 11 (9973-PD)

Images



Specifications

Property	Value
Material	Stainless Steel
Overall Length (mm)	125
To Fit Blades	11 & 15

Other Details

- This product can be viewed online at http://www.eyeline.org/products/view/bard_parker.

Punctum Dilator (9974D)

Images



Specifications

Property	Value
Material	Stainless Steel
Overall Length (mm)	101
Handle Style	Satin, Banded
Handle Diameter (mm)	4

Other Details

- This product can be viewed online at http://www.eyeline.org/products/view/punctum_dilator.

Lacrimal Probes (Set of 3) (9975D)

Images



Specifications

Property	Value
Material	Stainless Steel
Overall Length (mm)	120
Wire Diameter (mm)	1, 0.7, 0.5
Labelling	3/4, 1/2, 00/0

Other Details

- This product can be viewed online at http://www.eyeline.org/products/view/lacrimal_probes.