


Micro spatula (Tayyib) (9940D)

Images



Specifications

Property	Value
Material	Stainless Steel/Alloy
Overall Length (mm)	127
Tip Width (mm)	0.4
Tip Thickness (of Single Tip) (mm)	0.2
Handle Style	Knurled
Handle Diameter (mm)	5.5
Handle to First Bend Length (mm)	10
First Bend to Second Bend Length (mm)	11
Tip Style	Curved Micro-Spatula
Tip Length (mm)	1
Orientation (Left/Right/Either Hand)	Either Hand

Other Details


- This product can be viewed online at http://www.eyeline.org/products/view/microspatula_tayyib.

Micro Finger (9941D)

Images



tip detail



Specifications

Property	Value
Material	Stainless Steel/Alloy
Overall Length (mm)	127
Tip Width (mm)	0.4
Tip Thickness (of Single Tip) (mm)	0.4
Handle Style	Knurled
Handle Diameter (mm)	5.5
Handle to First Bend Length (mm)	10
First Bend to Second Bend Length (mm)	11
Tip Style	Straight Micro-Bar
Tip Length (mm)	2
Orientation (Left/Right/Either Hand)	Either Hand

Other Details


- This product can be viewed online at <http://www.eyeline.org/products/view/microfinger>.

Nagahara Chopper (9942D)

Images



tip detail



Specifications

Property	Value
Material	Stainless Steel/Alloy
Overall Length (mm)	127
Tip Width (mm)	0.4
Handle Style	Knurled
Handle Diameter (mm)	5.5
Handle to First Bend Length (mm)	10
First Bend to Second Bend Length (mm)	11
Tip Style	Blade Aligned with Shaft
Tip Length (mm)	1.25
Orientation (Left/Right/Either Hand)	Either Hand

Other Details


- This product can be viewed online at http://www.eyeline.org/products/view/chopper_nagahara.

Green's Chopper (for Right-Handed Surgeon) (9943D)

Images



tip detail



Specifications

Property	Value
Material	Stainless Steel/Alloy
Overall Length (mm)	127
Tip Width (mm)	0.3
Tip Thickness (of Single Tip) (mm)	0.3
Handle Style	Knurled
Handle Diameter (mm)	5.5
Handle to First Bend Length (mm)	10
First Bend to Second Bend Length (mm)	10
Tip Style	Green's Style Multi-Curved
Tip Length (mm)	2
Orientation (Left/Right/Either Hand)	Left Hand

Other Details


- This product can be viewed online at http://www.eyeline.org/products/view/chopper_green_right.

Green's Chopper (for Left-Handed Surgeon) (9944D)

Images



tip detail



Specifications

Property	Value
Material	Stainless Steel/Alloy
Overall Length (mm)	127
Tip Width (mm)	0.3
Tip Thickness (of Single Tip) (mm)	0.3
Handle Style	Knurled
Handle Diameter (mm)	5.5
Handle to First Bend Length (mm)	10
First Bend to Second Bend Length (mm)	10
Tip Style	Green's Style Multi-Curved
Tip Length (mm)	2
Orientation (Left/Right/Either Hand)	Right Hand

Other Details


- This product can be viewed online at http://www.eyeline.org/products/view/chopper_green_left.

Sinskey Implant Lens Dialler (9945D)

Images



tip detail



Description

Specifications

Property	Value
Material	Stainless Steel/Alloy
Overall Length (mm)	127
Tip Width (mm)	0.2
Tip Thickness (of Single Tip) (mm)	0.2
Handle Style	Knurled
Handle Diameter (mm)	5.5
Handle to First Bend Length (mm)	10
First Bend to Second Bend Length (mm)	11
Tip Style	Plain Dialler
Tip Length (mm)	0.9
Orientation (Left/Right/Either Hand)	Either Hand

Other Details


- This product can be viewed online at http://www.eyeline.org/products/view/sinskey_dialler.

Lester T-Dialler (9946D)

Images



tip detail



Specifications


Property	Value
Material	Stainless Steel/Alloy
Overall Length (mm)	127
Tip Width (mm)	0.40
Tip Thickness (of Single Tip) (mm)	0.4
Handle Style	Knurled
Handle Diameter (mm)	5.5
Handle to First Bend Length (mm)	10
First Bend to Second Bend Length (mm)	11
Tip Style	T-Shaped, Bulb on Tip
Tip Length (mm)	1.4
Orientation (Left/Right/Either Hand)	Either Hand

Other Details


- This product can be viewed online at http://www.eyeline.org/products/view/lester_dialler.

Bechert Y-Dialler (9947D)

Images



tip detail



Specifications


Property	Value
Material	Stainless Steel/Alloy
Overall Length (mm)	127
Tip Width (mm)	0.7
Handle Style	Knurled
Handle Diameter (mm)	5.5
Handle to First Bend Length (mm)	10
First Bend to Second Bend Length (mm)	11
Tip Style	Forked Tip, Rounded Ends
Tip Length (mm)	0.5
Orientation (Left/Right/Either Hand)	Either Hand

Other Details

- This product can be viewed online at http://www.eyeline.org/products/view/bechert_dialler.

Mushroom/Barrett Manipulator (9948D)

Images



Specifications



Property	Value
Material	Stainless Steel
Overall Length (mm)	127
Tip Width (mm)	0.65
Tip Thickness (of Single Tip) (mm)	0.65
Handle Style	Knurled
Handle Diameter (mm)	5.5
Handle to First Bend Length (mm)	10
First Bend to Second Bend Length (mm)	10
Tip Style	Mushroom
Tip Length (mm)	1.5
Orientation (Left/Right/Either Hand)	Either Hand

Other Details

- This product can be viewed online at http://www.eyeline.org/products/view/mushroom_barrett.

Iris Repositor flat blade (9949D)

Images



Specifications


Property	Value
Material	Stainless Steel/Alloy
Overall Length (mm)	127
Tip Width (mm)	0.2
Tip Thickness (of Single Tip) (mm)	0.5
Handle Style	Knurled
Handle Diameter (mm)	5.5
Handle to First Bend Length (mm)	10
First Bend to Second Bend Length (mm)	11
Tip Style	Angled Spatulate
Orientation (Left/Right/Either Hand)	Either Hand

Other Details

- This product can be viewed online at http://www.eyeline.org/products/view/irisrepositor_flat.

Micro paddle (Drysdale) (9939D)

Images



Specifications

Property	Value
Material	Stainless steel

Other Details

- This product can be viewed online at http://www.eyeline.org/products/view/micropaddle_drysdale.