


Bates Needle Holder (9935D)

Images



Specifications

Property	Value
Material	Stainless Steel
Overall Length (mm)	124
Tip - Shoulder Length (mm)	38
Tip + Handle Length (mm)	76
Handle Length (mm)	42
Heel - Handle (Spring) Length (mm)	48
Max Handle Width (mm)	7.9
Tip Width (mm)	0.7
Tip Thickness (of Single Tip) (mm)	0.4
Handle Style	Satin
Pivot - Handle (Shaft) Length (mm)	23
Tip - Pivot Length (mm)	15
Jaw Style	Long Curve
Effective Jaw Length (mm)	10
Shaft Width at Pivot (Blades Closed) (mm)	3.2
Shaft Thickness at Pivot (Combined) (mm)	2.6

Other Details

- This product can be viewed online at <http://www.eyeline.org/products/view/bates>.

Capsulorhexis Needle Bender (9937D)

Images



Specifications

Property	Value
Material	Stainless Steel
Overall Length (mm)	105
Tip Width (mm)	1.5
Tip Thickness (of Single Tip) (mm)	1
Tip - Pivot Length (mm)	18
Jaw Style	Plain Straight
Effective Jaw Length (mm)	15
Max Jaw Width (mm)	2.6

Other Details

- This product can be viewed online at http://www.eyeline.org/products/view/needle_bender.

Castroviejo Needle holder straight jaw (with lock) (9936D)

Images



Specifications

Property	Value
----------	-------

Other Details

- This product can be viewed online at http://www.eyeline.org/products/view/castroviejo_straight.

Castroviejo Needle holder curved jaw (with lock) (9938D)

Images



Specifications

Property	Value
----------	-------

Other Details

- This product can be viewed online at http://www.eyeline.org/products/view/castroviejo_curved.