

Conjunctival Scissors Curved Blunt (9928D)

Images



Specifications

Property	Value
Material	Stainless Steel
Overall Length (mm)	112
Pivot - Handle (Shaft) Length (mm)	51
Tip - Pivot Length (mm)	30
Blade Length (mm)	24
Effective Cut Length (mm)	24
Shaft Width at Pivot (Blades Closed) (mm)	5.8
Shaft Thickness at Pivot (Combined) (mm)	4.1
Radius of Blade Curve (mm)	100

Other Details

- This product can be viewed online at http://www.eyeline.org/products/view/conjunctival_curved_blunt.

Conjunctival Scissors Straight Pointed (9929D)

Images



Specifications

Property	Value
Material	Stainless Steel
Overall Length (mm)	112
Pivot - Handle (Shaft) Length (mm)	56
Tip - Pivot Length (mm)	30
Blade Length (mm)	24
Effective Cut Length (mm)	24
Shaft Width at Pivot (Blades Closed) (mm)	5.6
Shaft Thickness at Pivot (Combined) (mm)	2.8

Other Details

- This product can be viewed online at http://www.eyeline.org/products/view/conjunctival_straight_pointed.

Enucleation Scissors (9931D)

Images



Specifications

Property	Value
Material	Stainless Steel
Overall Length (mm)	114
Pivot - Handle (Shaft) Length (mm)	55
Tip - Pivot Length (mm)	31
Blade Length (mm)	28
Effective Cut Length (mm)	26
Shaft Width at Pivot (Blades Closed) (mm)	5.6
Shaft Thickness at Pivot (Combined) (mm)	3.6
Radius of Blade Curve (mm)	100

Other Details

- This product can be viewed online at <http://www.eyeline.org/products/view/enucleation>.

Dressing Scissors (9932D)

Images



Specifications

Property	Value
Material	Stainless Steel
Overall Length (mm)	130
Tip - Pivot Length (mm)	50
Blade Length (mm)	40
Effective Cut Length (mm)	40
Shaft Width at Pivot (Blades Closed) (mm)	9.8
Shaft Thickness at Pivot (Combined) (mm)	5.7
Radius of Blade Curve (mm)	200

Other Details

- This product can be viewed online at <http://www.eyeline.org/products/view/dressing>.