


# Bechert Y-Dialler (9947D)

## Images



tip detail



## Specifications

| Property                              | Value                    |
|---------------------------------------|--------------------------|
| Material                              | Stainless Steel/Alloy    |
| Overall Length (mm)                   | 127                      |
| Tip Width (mm)                        | 0.7                      |
| Handle Style                          | Knurled                  |
| Handle Diameter (mm)                  | 5.5                      |
| Handle to First Bend Length (mm)      | 10                       |
| First Bend to Second Bend Length (mm) | 11                       |
| Tip Style                             | Forked Tip, Rounded Ends |
| Tip Length (mm)                       | 0.5                      |
| Orientation (Left/Right/Either Hand)  | Either Hand              |

## Other Details

- This product can be viewed online at [http://www.eyeline.org/products/view/bechert\\_dialler](http://www.eyeline.org/products/view/bechert_dialler).