


# Green's Chopper (for Right-Handed Surgeon) (9943D)

## Images



tip detail



## Specifications

Property	Value
Material	Stainless Steel/Alloy
Overall Length (mm)	127
Tip Width (mm)	0.3
Tip Thickness (of Single Tip) (mm)	0.3
Handle Style	Knurled
Handle Diameter (mm)	5.5
Handle to First Bend Length (mm)	10
First Bend to Second Bend Length (mm)	10
Tip Style	Green's Style Multi-Curved
Tip Length (mm)	2
Orientation (Left/Right/Either Hand)	Left Hand

## Other Details

- This product can be viewed online at [http://www.eyeline.org/products/view/chopper\\_green\\_right](http://www.eyeline.org/products/view/chopper_green_right).