

Conjunctival Scissors Straight Pointed (9929D)

Images



Specifications

Property	Value
Material	Stainless Steel
Overall Length (mm)	112
Pivot - Handle (Shaft) Length (mm)	56
Tip - Pivot Length (mm)	30
Blade Length (mm)	24
Effective Cut Length (mm)	24
Shaft Width at Pivot (Blades Closed) (mm)	5.6
Shaft Thickness at Pivot (Combined) (mm)	2.8

Other Details

- This product can be viewed online at http://www.eyeline.org/products/view/conjunctival_straight_pointed.