

# Enucleation Scissors (9931D)

---

## Images



## Specifications

Property	Value
Material	Stainless Steel
Overall Length (mm)	114
Pivot - Handle (Shaft) Length (mm)	55
Tip - Pivot Length (mm)	31
Blade Length (mm)	28
Effective Cut Length (mm)	26
Shaft Width at Pivot (Blades Closed) (mm)	5.6
Shaft Thickness at Pivot (Combined) (mm)	3.6
Radius of Blade Curve (mm)	100

## Other Details

- This product can be viewed online at <http://www.eyeline.org/products/view/enucleation>.