

# Jayles Forceps (9918D)

---

## Images



## Specifications

Property	Value
Tooth Length + Tip Thickness (mm)	0.7
Tooth Configuration	1 into 2 protruding
Tooth Length (mm)	0.2
Tip Thickness (of Single Tip) (mm)	0.5
Tip Width (mm)	0.7
Max Handle Thickness (mm)	2
Max Handle Width (mm)	7.5
Heel - Handle (Spring) Length (mm)	32
Tip + Handle Length (mm)	80
Platform Length (mm)	7
Overall Length (mm)	112
Material	Stainless Steel

## Other Details

- This product can be viewed online at <http://www.eyeline.org/products/view/jayles>.