

McPherson Angled Forceps (9908D)

Images



Specifications

Property	Value
Material	Stainless Steel
Overall Length (mm)	110
Platform Length (mm)	7
Tip - Shoulder Length (mm)	25
Tip + Handle Length (mm)	78
Handle Length (mm)	53
Heel - Handle (Spring) Length (mm)	32
Max Handle Width (mm)	10
Max Handle Thickness (mm)	1.7
Tip Width (mm)	0.4
Tip Thickness (of Single Tip) (mm)	0.35
Number of Holes	3
Handle Style	Castroviejo
Handle Finish	Satin
Tip - Angle Length (mm)	8

Other Details

- This product can be viewed online at <http://www.eyeline.org/products/view/mcpherson>.