

# Capsulorhexis Needle Bender (9937D)

---

## Images



## Specifications

Property	Value
Material	Stainless Steel
Overall Length (mm)	105
Tip Width (mm)	1.5
Tip Thickness (of Single Tip) (mm)	1
Tip - Pivot Length (mm)	18
Jaw Style	Plain Straight
Effective Jaw Length (mm)	15
Max Jaw Width (mm)	2.6

## Other Details

- This product can be viewed online at [http://www.eyeline.org/products/view/needle\\_bender](http://www.eyeline.org/products/view/needle_bender).