

Corneal Plain Forceps (9900D)

Images



Specifications

Property	Value
Handle Finish	Satin
Number of Holes	3
Tip Thickness (of Single Tip) (mm)	0.35
Tip Width (mm)	0.4
Max Handle Thickness (mm)	1.7
Max Handle Width (mm)	10
Heel - Handle (Spring) Length (mm)	31
Handle Length (mm)	53
Tip + Handle Length (mm)	78
Tip - Shoulder Length (mm)	25.0
Overall Length (mm)	110
Material	Stainless Steel

Other Details

- This product can be viewed online at http://www.eyeline.org/products/view/plain_hoskins.