



Sinskey Implant Lens Dialler (9945D)

Images



tip detail



Description

Specifications

Property	Value
Material	Stainless Steel/Alloy
Overall Length (mm)	127
Tip Width (mm)	0.2
Tip Thickness (of Single Tip) (mm)	0.2
Handle Style	Knurled
Handle Diameter (mm)	5.5
Handle to First Bend Length (mm)	10
First Bend to Second Bend Length (mm)	11
Tip Style	Plain Dialler
Tip Length (mm)	0.9
Orientation (Left/Right/Either Hand)	Either Hand

Other Details

- This product can be viewed online at http://www.eyeline.org/products/view/sinskey_dialler.