

St Martins Forceps (9917D)

Images



Specifications

Property	Value
Tooth Length + Tip Thickness (mm)	0.7
Tooth Configuration	1 into 2
Tooth Length (mm)	0.27
Tip Thickness (of Single Tip) (mm)	0.5
Tip Width (mm)	0.77
Max Handle Thickness (mm)	2
Max Handle Width (mm)	11
Heel - Handle (Spring) Length (mm)	30
Handle Length (mm)	40
Tip + Handle Length (mm)	62
Tip - Shoulder Length (mm)	12
Platform Length (mm)	6
Overall Length (mm)	92
Material	Stainless Steel

Other Details

- This product can be viewed online at http://www.eyeline.org/products/view/st_martin.