

# Corneal Toothed Forceps with Platform (9906D)

---

## Images



## Specifications

Property	Value
Material	Stainless Steel
Overall Length (mm)	110
Platform Length (mm)	6
Tip - Shoulder Length (mm)	25
Tip + Handle Length (mm)	78
Handle Length (mm)	53
Heel - Handle (Spring) Length (mm)	32
Max Handle Width (mm)	10
Max Handle Thickness (mm)	1.7
Tip Width (mm)	0.3
Tip Thickness (of Single Tip) (mm)	0.35
Number of Holes	3
Handle Style	Castroviejo
Handle Finish	Satin
Tooth Length (mm)	0.12
Tooth Configuration	1 into 2

## Other Details

- This product can be viewed online at [http://www.eyeline.org/products/view/toothed\\_platforms](http://www.eyeline.org/products/view/toothed_platforms).