


# Corneal Tying Straight Forceps (9902D)

---

## Images



## Specifications

Property	Value
Material	Stainless Steel
Overall Length (mm)	110
Platform Length (mm)	6.0
Tip - Shoulder Length (mm)	20
Tip + Handle Length (mm)	80
Handle Length (mm)	60
Heel - Handle (Spring) Length (mm)	30
Max Handle Width (mm)	10
Max Handle Thickness (mm)	1.7
Tip Width (mm)	0.4
Tip Thickness (of Single Tip) (mm)	0.35
Number of Holes	3
Handle Style	Castroviejo
Handle Finish	Satin

## Other Details

- This product can be viewed online at [http://www.eyeline.org/products/view/tying\\_straight\\_long](http://www.eyeline.org/products/view/tying_straight_long).