

Utrata Curved Forceps (9910D)

Images



Specifications

Property	Value
Tip - Angle Length (mm)	13.0
Tooth Configuration	1 onto 1
Tooth Length (mm)	0.25
Handle Finish	Satin
Handle Style	Castroviejo
Number of Holes	3
Tip Thickness (of Single Tip) (mm)	0.25
Max Handle Thickness (mm)	1.7
Max Handle Width (mm)	10
Heel - Handle (Spring) Length (mm)	30
Tip + Handle Length (mm)	80
Overall Length (mm)	110
Material	Stainless Steel

Other Details

- This product can be viewed online at http://www.eyeline.org/products/view/ultrata_curved.