

Vannas Curved Scissors (9920D)

Images



Specifications

Property	Value
Material	Stainless Steel
Overall Length (mm)	82
Tip + Handle Length (mm)	33
Heel - Handle (Spring) Length (mm)	30
Max Handle Width (mm)	7.6
Max Handle Thickness (mm)	1.5
Tip Width (mm)	0.35
Tip Thickness (of Single Tip) (mm)	0.30
Handle Style	Ridged
Pivot - Handle (Shaft) Length (mm)	23
Tip Thickness (of Both Tips, Closed) (mm)	0.43
Tip - Pivot Length (mm)	10.4
Blade Length (mm)	8
Effective Cut Length (mm)	5.2
Shaft Width at Pivot (Blades Closed) (mm)	1.83
Shaft Thickness at Pivot (Combined) (mm)	1.72
Radius of Blade Curve (mm)	25

Other Details

- This product can be viewed online at http://www.eyeline.org/products/view/vannas_curved.