

# Vannas Straight Scissors (9921D)

## Images



## Specifications

Property	Value
Material	Stainless Steel
Overall Length (mm)	90
Tip + Handle Length (mm)	53
Heel - Handle (Spring) Length (mm)	37
Max Handle Width (mm)	7.7
Max Handle Thickness (mm)	1.5
Tip Width (mm)	0.8
Tip Thickness (of Single Tip) (mm)	0.25
Handle Style	Ridged
Pivot - Handle (Shaft) Length (mm)	24
Tip Thickness (of Both Tips, Closed) (mm)	0.57
Tip - Pivot Length (mm)	11
Blade Length (mm)	9.5
Effective Cut Length (mm)	6.5
Shaft Width at Pivot (Blades Closed) (mm)	1.9
Shaft Thickness at Pivot (Combined) (mm)	1.88

## Other Details

- This product can be viewed online at [http://www.eyeline.org/products/view/vannas\\_straight](http://www.eyeline.org/products/view/vannas_straight).