


Williamson Noble Scissors Curved Blunt (9925D)

Images



Specifications

Property	Value
Material	Stainless Steel
Overall Length (mm)	123
Tip + Handle Length (mm)	80
Heel - Handle (Spring) Length (mm)	43
Max Handle Width (mm)	7.9
Max Handle Thickness (mm)	2.4
Tip Width (mm)	0.66
Tip Thickness (of Single Tip) (mm)	0.40
Handle Style	Ridged
Pivot - Handle (Shaft) Length (mm)	25
Tip Thickness (of Both Tips, Closed) (mm)	0.95
Tip - Pivot Length (mm)	24
Blade Length (mm)	18
Effective Cut Length (mm)	14.8
Shaft Width at Pivot (Blades Closed) (mm)	4
Shaft Thickness at Pivot (Combined) (mm)	3.5
Radius of Blade Curve (mm)	250

Other Details

- This product can be viewed online at http://www.eyeline.org/products/view/williamson_noble_curved_blunt.