


Williamson Noble Scissors Curved Pointed (9924D)

Images



Specifications

| Property | Value |
|---|-----------------|
| Material | Stainless Steel |
| Overall Length (mm) | 125 |
| Tip + Handle Length (mm) | 82 |
| Heel - Handle (Spring) Length (mm) | 43 |
| Max Handle Width (mm) | 7.9 |
| Max Handle Thickness (mm) | 2.4 |
| Tip Width (mm) | 0.66 |
| Tip Thickness (of Single Tip) (mm) | 0.40 |
| Handle Style | Ridged |
| Pivot - Handle (Shaft) Length (mm) | 25 |
| Tip Thickness (of Both Tips, Closed) (mm) | 0.95 |
| Tip - Pivot Length (mm) | 24.2 |
| Blade Length (mm) | 18.2 |
| Effective Cut Length (mm) | 15 |
| Shaft Width at Pivot (Blades Closed) (mm) | 4 |
| Shaft Thickness at Pivot (Combined) (mm) | 3.5 |
| Radius of Blade Curve (mm) | 250 |

Other Details

- This product can be viewed online at http://www.eyeline.org/products/view/williamson_noble_curved_pointed

